## **Cleaning and Sterilization Instructions**

The following recommendations are for processing Stryker reusable medical devices. Reusable devices include the nylon carving block, silicone and stainless steel sizing sets and fixation kits (including titanium plates and screws, surgical instruments, instrument trays and lids). **THESE RECOMMENDATIONS DO NOT APPLY TO MEDPOR POLYETHYLENE IMPLANTS.** MEDPOR polyethylene implants are provided Sterile for Single Use Only and should not be re-sterilized under any conditions.

## Manual Cleaning Instructions

1. Pre-wash in warm tap water for 2 minutes.
2. Completely submerge the reusable devices in an enzyme solution and allow to soak for 20 minutes. Use a soft-bristled, nylon brush to gently scrub the device until all visible soil has been removed. Particular attention must be given to crevices, lumens, mated surfaces, connected and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft-bristled brush (i.e. pipe cleaner brush).
3. Remove the device from the enzyme solution and rinse in warm tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, blind holes and other difficult-to-reach areas.
4. Place prepared cleaning agents in a sonication unit. Completely submerge devices in cleaning solution and sonicate for 10 minutes at 45 – 50 kHz.
5. Rinse instruments in Sterile water for at least 3 minutes or until there is no sign of blood or soil on the device or in the rinse stream. Thoroughly and aggressively flush lumens, blind holes and other difficult-to-reach areas. Repeat sonication and rinse steps above.
6. Remove excess moisture and clean with an absorbent and non-shedding wipe.

*Note: If stainless steel instruments are stained or corroded, an acidic, anti-corrosion agent in ultrasonic cleaner may be sufficient to remove surface deposits. Care must be taken to thoroughly rinse acid from devices. Acidic, anti-corrosion agents should only be used on an as-needed basis.*

## Combination Manual / Automated Cleaning Instructions

1. Pre-wash in warm tap water for 2 minutes.
2. Completely submerge the reusable devices in an enzyme solution and allow to soak for 20 minutes. Use a soft-bristled, nylon brush to gently scrub the device until all visible soil has been removed. Particular attention must be given to crevices, lumens, mated surfaces, connected and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft-bristled brush (i.e. pipe cleaner brush).
3. Remove the device from the enzyme solution and rinse in warm tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, blind holes and other difficult-to-reach areas.
4. Place prepared cleaning agents in a sonication unit. Completely submerge devices in cleaning solution and sonicate for 10 minutes at 45 – 50 kHz.
5. Rinse instruments in sterile water for at least 3 minutes or until there is no sign of blood or soil on the device or in the rinse stream. Thoroughly and aggressively flush lumens, blind holes and other difficult-to-reach areas. Repeat sonication and rinse steps above.
6. Place device in a suitable washer / disinfectant basket and process through a standard instrument cycle. The following minimum parameters are essential for thorough cleaning and disinfection:

Step	Description
1.	2 minutes pre-wash with warm tap water
2.	20 seconds enzyme spray with hot tap water
3.	1 minute enzyme soak
4.	15 seconds cold tap water rinse (X2)
5.	2 minutes detergent wash with hot tap water (64-66°C/146-150°F)
6.	15 seconds hot tap water rinse
7.	2 minutes thermal rinse (80-93°C/176-200°F)
8.	10 seconds purified water rinse with optional lubricant (64-66°C/146-150°F)
9.	7-30 minutes hot air dry (>=116°C/240°F)

*Note: The washer / disinfectant manufacturer's instruction should be strictly adhered to.*

## STERILE PACKAGING Packaging Individual Instruments, Kits or Silicone Devices

1. Commercially available, medical grade, steam sterilization wraps or pouches (paper, Tyvek or equivalent) may be used to double package single instruments. Ensure the package is large enough to contain the instrument without stressing the seals.
2. Only intact components or instruments should be processed. Functionality must be verified following cleaning and prior to packaging.

## STERILIZATION

1. Do not re-sterilize MEDPOR polyethylene implants or any other Single Use Only components.
2. These sterilization instructions are consistent with AORN and ANSI / AAMI / ISO guidelines. They should be used for items supplied non sterile (excluding non-sterile sample products provided for demonstration purposes only), for reprocessing reusable devices, or for sterile items that can be re-sterilized.
3. Where possible, reusable instruments should be disassembled for sterilization and all surfaces must be accessible to the sterilizing medium.
4. Instruments, silicone sizer sets, carving blocks, fixation kits including titanium plates and screws, surgical instruments, instrument trays and lids may be wrapped in a standard medical grade, steam sterilization wrap using the AAMI double wrap method or equivalent.
5. Do not sterilize / re-sterilize devices in the protective bag or container in which they are supplied.
6. Precaution/Warning: Components should be sterilized loose and not stacked in layers. Follow the sterilizer manufacturer instructions for loading patterns.

All reusable instruments, instrument kits, fixation kits and carving blocks should be steam sterilized using these parameters:

Type	Minimum Temperature	Minimum Exposure <sup>6</sup>	Minimum Dry Time <sup>7</sup>
U.K. Pre-Vacuum/ Pulsating Vacuum <sup>1,3,4</sup>	134°C / 273°F	15 minutes	30 minutes
U.S. Pre-Vacuum/ Pulsating Vacuum <sup>2,3,5</sup>	132°C / 270°F	15 minutes	Varies by load configuration and sterilizer type
Gravity Cycle	121°C / 254 °F	60 minutes	

All reusable silicone sizer sets should be steam sterilized using these parameters:

Type	Minimum Temperature	Minimum Exposure <sup>6</sup>	Minimum Dry Time <sup>7</sup>
U.K. Pre-vacuum/ Pulsating Vacuum <sup>1,3,4</sup>	134°C / 273°F	15 minutes	30 minutes
U.S. Pre-Vacuum/ Pulsating Vacuum <sup>2,3</sup>	132°C / 270°F	15 minutes	Varies by load configuration and sterilizer type
Gravity Cycle	121°C / 254 °F	30 minutes	

1 Minimum validated steam sterilization time required to achieve a 10<sup>-6</sup> sterility assurance level (SAL).

2 Minimum validated steam sterilization temperature required to achieve a 10<sup>-6</sup> sterility assurance level (SAL).

3 Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in this table.

4 Disinfection / steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CJD contamination.

5 For universal instrument cases without defined load configurations.

6 AAMI / AORN steam sterilization cycles with longer cycle times than those listed are also acceptable.

7 Drying times vary according to load size and should be increased for larger loads.

Note: Instructions above are based on validation studies. The recommended maximum number of re-sterilization cycles endorsed by any one item should be no greater than five.

Copyright © 2010 Stryker. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or services: **MEDPOR** and **Stryker**. All other trademarks are trademarks of their respective owners or holders. ECO-069-97A

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Copyright © 2010 Stryker. MEDPOR es una marca comercial registrada de Stryker Corporation. La marca de servicio Stryker es propiedad de Stryker Corporation o se utiliza bajo su licencia y autorización. ECO-069-97A
Copyright © 2010 Stryker. MEDPOR ist ein eingetragenes Warenzeichen der Stryker Corporation. Die Stryker Servicemarke ist entweder Eigentum der Stryker Corporation oder wird unter ihrer Lizenz und Aufsicht verwendet. ECO-069-97A
Copyright © 2010 Stryker. MEDPOR est une marque déposée de Stryker Corporation. La marque de service Stryker appartient à, ou est utilisée sous licence et autorité de, Stryker Corporation. ECO-069-97A