
## COI™ ILLUSTRATIVE GUIDE *MEDPOR® Surgical Implant*

ILLUSTRATIVE GUIDE  
Product Insert Addendum for Conical Orbital Implant

The successful utilization of MEDPOR Conical Orbital Implants is technique-dependent.

## Labeling Symbols

Attention, See Instructions for Use

Do Not Reuse

Sterilized using Ethylene Oxide

Sterilized using Irradiation

Use by Date

Legal Manufacturer

Authorized Representative in the European Community

Batch Code / Lot Number

Catalog Number

Keep Dry

Not Sterile


Do Not Use if Package is Damaged

## Conical Orbital Implant Designed in Conjunction with Peter A.D. Rubin, M.D.

The Conical Orbital Implant (COI) was designed and developed to address many of the common problems associated with the correction of the anophthalmic socket.

Unique design elements have been incorporated into the overall conical shape, including a superior projection and channels for the rectus muscles. Three sizes are available. Catalog number 7222, Conical Orbital Implant - 4ml, Catalog number 7223, Conical Orbital Implant - 5ml and Catalog number 80026, Conical Orbital Implant - 6ml. The redistributed volume of the superior aspect of the implant helps minimize the post operative superior sulcus defect.

This shape is intended for use after standard enucleation procedures. The implant is primarily intended to (1) fill the void volume resulting from an enucleated eye, (2) provide a method for reattaching the rectus muscles, and (3) provide a compatible surface for an overlying ocular prosthesis, with or without the use of the Motility Coupling Post (MCP).


## Features


MEDPOR® Implants are composed of biocompatible porous polyethylene for use in orbital reconstruction procedures. The interconnecting, omni-directional pore structure of the implant allows for rapid vascularization and tissue ingrowth.

The superior projection redistributes volume to the upper portion of the implant to minimize post operative superior sulcus defects.

The relative flat anterior surface provides adequate depth for the precision fitting of an ocular prosthesis.

Muscle channels provide location and access for attachment of the lateral, medial, inferior and superior rectus muscles. The superior channel has preformed suture holes to facilitate rectus muscle placement.



**Fig. 1** In smaller orbits conjunctival shortage, implant modification with a scalpel can be performed. The extraocular muscle channels permit central and anterior placement of the extraocular muscles.

**Fig. 2** A fascial or scleral graft may be sutured to the anterior face.

**Fig. 7** Straighten a needle in order to pass it through the suture holes for the superior rectus muscle.

**Fig. 8** Slide the implant down into the socket in a "cable car" fashion utilizing the superior rectus muscle sutures.


**Fig. 3** The anterior face of the implant is shown with an autologous fascial graft sutured in position.


**Fig. 4** A posterior view shows the unwrapped portion of the implant as well as the muscle channels.


**Fig. 9** Once placed, digital pressure will position the implant posteriorly in the orbit. Tie off the superior muscle sutures.


**Fig. 10** Bring the remaining muscle sutures under and through the fascia covering and tie them off about 5mm from the center of the anterior surface.


**Fig. 5** Following enucleation, posterior Tenon's is opened widely to permit deeper positioning of the implant in the orbit.


**Fig. 6** Prior to implantation the conjunctival edges are grasped and retracted to ease insertion.

**Fig. 11** The conjunctival edges are grasped and pulled anteriorly to facilitate closure.

**Fig. 12** Anterior Tenon's is closed meticulously without tension using interrupted sutures (i.e. 5-0 polyglactin [Vicryl™]).

**Fig. 13** The conjunctiva is closed without tension with a running suture (i.e. fast-absorbing plain-gut 6-0).

**Fig. 14** A conformer is positioned within the conjunctival cul-de-sac. A suture tarsorrhaphy may be beneficial.

**Fig. 15** If a tarsorrhaphy is placed, only a light pressure patch is needed.

**Fig. 16** Postoperative appearance of enucleated patient.

**Fig. 17** Patient fitted with ocular prosthesis.

## References

- Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Choi, J.C., Iwamoto, M.A., Bstandig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. *Ophthalmology* 101:1688-1696, 1994.
- Rubin, P.A.D. Conical porous polyethylene enucleation implant. *Journal of the American Society of Ocularists*, 25th Edition, 1994.
- Rubin, P.A.D., Blysk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

## Ordering Information

Implants are provided STERILE, packaged individually in double peel pouches.

Catalog No.	Description
7222	Conical Orbital Implant - 4ml
7223	Conical Orbital Implant - 5ml
80026	Conical Orbital Implant - 6ml

Both the 4ml and 5ml models have 18mm frontal dimensions with the 5ml having a 24 mm body length vs. 22mm in the 4ml model. The 6ml model has a 19mm frontal dimension with a body length 25mm. Surgeons should utilize proper surgical techniques and their clinical experience to determine appropriate procedures. Successful implantations are technique sensitive. Sound surgical judgment should be used in the selection, shaping, handling, and implantation of all MEDPOR® shapes. Please be familiar with the MEDPOR® Surgical Implant Product Information insert.

Call to order: 1-800-962-6558

U.S. Patent No. 5,466,373

Copyright © 2011, Howmedica Osteonics Corp.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or services: COI, MEDPOR, Stryker. All other trademarks are trademarks of their respective owners or holders.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Approximate	Approx	Lateral Projection	L.P.
Case	Cs	Length	Lnth
Centimeter	cm	Milliliter	ml
Contoured	Cont'd	Millimeter	mm
Cubic Centimeters	cc	Pack	pk
Diameter	DIA	Package	Pkg
Dimensions	DIM	Quantity	Qty
External	Ext	Thickness	Thick
French	Fr	Volume	vol
Inch	in	With	W/
Including	incl	Without	W/O
Inferior	Infer		

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## Bibliographie

- Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Choi, J.C., Iwamoto, M.A., Bständig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. Ophthalmology 101:1688-1696, 1994.
- Rubin, P.A.D. Conical porous polyethylene enucleation implant. Journal of the American Society of Ocularists, 25th Edition, 1994.
- Rubin, P.A.D., Blyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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**Abb. 17** Patient mit Augenprothese.

## Literatur

- Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Choi, J.C., Iwamoto, M.A., Bstendig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. *Ophthalmology* 101:1688-1696, 1994.
- Rubin, P.A.D. Conical porous polyethylene enucleation implant. *Journal of the American Society of Ocularists*, 25th Edition, 1994.
- Rubin, P.A.D., Bilyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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## Bibliografia

- Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Choi, J.C., Iwamoto, M.A., Bstandig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. *Ophthalmology* 101:1688-1696, 1994.
- Rubin, P.A.D. Conical porous polyethylene enucleation implant. *Journal of the American Society of Ocularists*, 25th Edition, 1994.
- Rubin, P.A.D., Blyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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- Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Choi, J.C., Iwamoto, M.A., Bständig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. *Ophthalmology* 101:1688-1696, 1994.
- Rubin, P.A.D. Conical porous polyethylene enucleation implant. *Journal of the American Society of Ocularists*, 25th Edition, 1994.
- Rubin, P.A.D., Blyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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## Referenser

- Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Choi, J.C., Iwamoto, M.A., Bstandig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. *Ophthalmology* 101:1688-1696, 1994.
- Rubin, P.A.D. Conical porous polyethylene enucleation implant. *Journal of the American Society of Ocularists*, 25th Edition, 1994.
- Rubin, P.A.D., Bilyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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Choi, J.C., Iwamoto, M.A., Bstendig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.

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Rubin, P.A.D., Blyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic

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**Figur 17** Patient udstyret med okulær protese.

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- Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Choi, J.C., Iwamoto, M.A., Bstendig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. Ophthalmology 101:1688-1696, 1994.
- Rubin, P.A.D. Conical porous polyethylene enucleation implant. Journal of the American Society of Ocularists, 25th Edition, 1994.
- Rubin, P.A.D., Bilyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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Choi, J.C., Iwamoto, M.A., Bstangid, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.

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Rubin, P.A.D., Blyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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## Referanslar

Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.

Choi, J.C., Iwamoto, M.A., Bstandig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.

Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. Ophthalmology 101:1688-1696, 1994.

Rubin, P.A.D. Conical porous polyethylene enucleation implant. Journal of the American Society of Ocularists, 25th Edition, 1994.

Rubin, P.A.D., Bilyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic

Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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## Bibliografia

- Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Choi, J.C., Iwamoto, M.A., Bstandig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. Ophthalmology 101:1688-1696, 1994.
- Rubin, P.A.D. Conical porous polyethylene enucleation implant. Journal of the American Society of Ocularists, 25th Edition, 1994.
- Rubin, P.A.D., Bilyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium; American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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- Choi, J.C., Rubin, P.A.D., Popham, J.K., Iwamoto, M.A., Shore, J.W. Presentation - MEDPOR® Motility Coupling Post: A New Coupling Device for MEDPOR® Orbital Implants. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Choi, J.C., Iwamoto, M.A., Bstendig, S., Rubin, P.A.D., Shore, J.W. Presentation - ASOPRS Thesis: MEDPOR® Motility Coupling Post: A Rabbit Model. American Society of Oculoplastic and Reconstructive Surgeons Annual Meeting, Chicago, IL, October, 1996.
- Karesh, J.W. and Dresner, S.C. High density porous polyethylene (MEDPOR) as a successful anophthalmic socket implant. *Ophthalmology* 101:1688-1696, 1994.
- Rubin, P.A.D. Conical porous polyethylene enucleation implant. *Journal of the American Society of Ocularists*, 25th Edition, 1994.
- Rubin, P.A.D., Bilyk, J.R., Popham, J.K., Shore, J.W. A New Enucleation Implant: Conical Shape with Sulcus Support. 25th Annual Scientific Symposium: American Society of Ophthalmic Plastic and Reconstructive Surgery, Insights & Innovations, October, 1994, pp. 43.

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A subsidiary of Stryker Corporation